



# R10 Glovebox Technique

*Jed Januch & Keven McDermott*

**Asbestos Site Evaluation, Communication and Cleanup**

**Keystone, Colorado  
September 22—26, 2003**



# What is it?

- Glovebox technique used to determine if a matrix that contains asbestos will release fibers when that matrix is disturbed
- Glovebox technique first used by R10 during vermiculite garden products study in 2000



# Why the Glovebox?

- R10 needed a quick and inexpensive method to determine whether asbestos in vermiculite or soil would become airborne when disturbed





# What is the Glovebox?



# Glovebox Specs

- Kewaunee Scientific Equipment
  - 4ft long, 3 ft high, 21 - 27 inches deep
  - volume approximately 680 liters
  - stainless steel w/ glass front
  - three portals to insert gloves
  - one large opening on side



# Air Monitoring Equipment Used in Glovebox

- 2 Allegro A100 high volume sampling pumps
  - air monitoring cassettes connected to pumps with tygon tubing
  - make-up air supplied through tygon tubing





# What did it look like?



# What was the process?

- Soil samples collected from asbestos contaminated site
- Some samples dried
- Samples placed in glovebox and sieved
- Samples agitated until dust was visible
- Air samples collected after agitation






# More Process

- Optimum loading of filter cassettes determined by PCM
  - duration and flow of air sampling
    - 12 minutes @ 2.5 - 3.5 liters/min
  - waiting time after agitation of samples
    - one hour







# More Process

- Duplicate set of samples collected under optimum flow and duration and analyzed using TEM





# Advantages of Glovebox

- The price was right - surplus equipment
- Airtight - potential exposure to asbestos minimized
- Easy to equip with tubing and pumps



# Disadvantages of Glovebox

- Difficult to clean
  - HEPA vac after each sample
  - rinse with deionized water
  - dry with fiber free wipes
  - time consuming







# Why We Like the R10 Glovebox

- Tells us quickly whether fibers present in samples become airborne



# Glovebox Limitations

- Artificial conditions not necessarily representative of “real world” exposures

