

APPENDIX B

SUPPLEMENTAL LIST OF ACUTE TOXICITY TEST SPECIES

TEST ORGANISM		TEST TEMP (°C)	LIFE STAGE
FRESHWATER SPECIES: VERTEBRATES - WARMWATER			
<i>Cyprinella leedsi</i> ¹	Bannerfin shiner	25	1-14 days
<i>Lepomis macrochirus</i>	Bluegill sunfish	20,25	" "
<i>Ictalurus punctatus</i>	Channel catfish	"	" "
FRESHWATER SPECIES: INVERTEBRATES - COLDWATER			
<i>Pteronarcys spp.</i>	Stoneflies*	12	larvae
<i>Pacifastacus leniusculus</i>	Crayfish*	"	juveniles
<i>Baetis spp.</i>	Mayflies*	"	nymphs
<i>Ephemerella spp.</i>	"	"	"
FRESHWATER SPECIES: INVERTEBRATES - WARMWATER			
<i>Hyaella spp.</i>	Amphipods	20,25	juveniles
<i>Gammarus lacustris</i>	"	"	"
<i>G. fasciatus</i>	"	"	"
<i>G. pseudolimnaeus</i>	"	"	"
<i>Hexagenia limbata</i>	Mayflies	"	nymphs
<i>H. bilineata</i>	"	"	"
<i>Chironomus spp.</i>	Midges	"	larvae

*Stoneflies, crayfish, and mayflies may have to be field collected and acclimated for a period of time to ensure the health of the organisms and that stress from collection is past. Species identification must be verified.

¹ Test conditions for *Cyprinella leedsi* and *Holmesimysis costata* are found in Table 14 and Table 19, respectively, in Section 9.

SUPPLEMENTAL LIST OF ACUTE TOXICITY TEST SPECIES (CONTINUED)

TEST ORGANISM		TEST TEMP (°C)	SALINITY (‰)	LIFE STAGE
MARINE AND ESTUARINE SPECIES: VERTEBRATES - COLDWATER				
<i>Parophrys vetulus</i>	English sole	12	32-34	1-90 days
<i>Citharichys sitigmaeus</i>	Sanddab	"	"	" "
<i>Pseudopleuronectes americanus</i>	Winter flounder	"	"	post metamorphosis
MARINE AND ESTUARINE SPECIES: VERTEBRATES - WARMWATER				
<i>Paralichthys dentatus</i>	Flounder	20,25	32-34	1-90 days
<i>P. lethostigma</i>	"	"	"	" "
<i>Fundulus simillis</i>	Killifish	"	20-32	1-30 days
<i>Fundulus heteroclitus</i>	Mummichog	"	25-32	" "
<i>Lagodon rhomboides</i>	Pinfish	"	20-32	1-90 days
<i>Orthipristis chrysoptera</i>	Pigfish	"	15-30	" "
<i>Leostomus xanthurus</i>	Spot	"	10-30	" "
<i>Gasterosteus aculeatus</i>	Threespine stickleback	"	20-32	1-30 days
<i>Atherinops affinis</i>	Topsmelt	21	10-30	7-15 days
MARINE AND ESTUARINE SPECIES: INVERTEBRATES - COLDWATER				
<i>Pandalus jordani</i>	Oceanic shrimp	12	25-32	juvenile
<i>Strongylocentrotus droebachiensis</i>	Green sea urchin	"	32-34	gametes/embryo
<i>Strongylocentrotus purpuratus</i>	Purple sea urchin	"	"	" "
<i>Dendraster excentricus</i>	Sand dollar	"	"	" "
<i>Cancer magister</i>	Dungeness crab	"	"	juvenile
<i>Holmesimysis costata</i> ²	Mysid	"	"	1-5 days
MARINE AND ESTUARINE SPECIES: INVERTEBRATES - WARMWATER				
<i>Callinectes sapidus</i>	Blue crab	20,25	10-30	juvenile
<i>Palaemonetes pugio</i>	Grass shrimp	"	10-32	1-10 days
<i>P. vulgaris</i>	" "	"	"	" "
<i>P. intermedius</i>	" "	"	"	" "
<i>Penaeus setiferus</i>	White shrimp	"	20-32	post-larval
<i>Penaeus duorarum</i>	Pink shrimp	"	"	" "
<i>Penaeus aztecus</i>	Brown shrimp	"	"	" "
<i>Crangon septemspinosa</i>	Sand shrimp	"	25-32	" "
<i>Mysidopsis almyra</i>	Mysid	"	10-32	1-5 days
<i>Neomysis americana</i>	"	"	"	" "
<i>Metamysidopsis elongata</i>	"	"	"	" "
<i>Crassostrea virginica</i>	American oyster	"	20-32	embryo
<i>Crassostrea gigas</i>	Pacific oyster	"	25-32	"
<i>Arbacia punctulata</i>	Purple sea urchin	"	32-34	gametes/embryo

² Test conditions for *Holmesimysis costata* are found in Table 19.

APPENDIX C

DILUTOR SYSTEMS

Two proportional dilutor systems are illustrated: the solenoid valve system, and the vacuum siphon system.

1. Solenoid and Vacuum Siphon Dilutor Systems

The designs of the solenoid and vacuum siphon dilutor systems incorporate features from devices developed by many other Federal and state programs, and have been shown to be very versatile for on-site bioassays in mobile laboratories, as well as in fixed (central) laboratories. The Solenoid Valve system is fully controlled by solenoids (Figures 1, 2, and 3), and is preferred over the vacuum siphon system. The Vacuum Siphon system (Figures 1, 4, and 5), however, is acceptable. The dilution water, effluent, and pre-mixing chambers for both systems are illustrated in Figures 6, 7, and 8. Both systems employ the same control panel (Figure 9).

If in the range-finding test, the LC50 of the effluent falls in the concentration range, 6.25-100%, pre-mixing is not required. The pre-mixing chamber is bypassed by running a TYGON® tube directly from the effluent in-flow pipe to chamber E-2 (see Figures 3 and 5), and Chambers E-1 and D-1 and the pre-mixing chamber are deactivated.

The dilutor systems described here can also be used to conduct tests of the toxicity of pure compounds by equipping the control panel with an auxiliary power receptacle to operate a metering pump to deliver an aliquot of the stock solution of the pure compound directly to the mixing chamber during each cycle. In this case, chamber E-1 is de-activated and chamber D-1 is calibrated to deliver a volume of 2000 mL, which is used to dilute the aliquot to the highest concentration used in the toxicity test.




Figure 1. Photographs of the solenoid valve system (left), and the vacuum siphon system (right).





Figure 2. Solenoid valve dilutor system, general diagram (not to scale).



NOTE: WHEN 100% EFFLUENT IS USED AS THE HIGHEST EFFLUENT CONCENTRATION, E-1, D-1, AND THE PRE-MIXING CHAMBER ARE BYPASSED BY CONNECTING A TYGON TUBE TO THE EFFLUENT INFLOW, AND RUNNING IT DIRECTLY TO E-2. IN THIS CASE, SOLENOIDS FOR E-1 AND D-1, AND THE PRE-MIXING CHAMBER ARE DISCONNECTED. D-2 + E-3 = 80% EFFLUENT; D-3 + E-4 = 85% EFFLUENT, ETC.

Figure 3. Solenoid valve dilutor system, detailed diagram (not to scale).

SOLENOID SYSTEM EQUIPMENT LIST

1. Dilator Glass.
2. Stainless Steel Solenoid Valves
 - a. 3, normally open, two-way, 55 psi, water, 1/4" pipe size, 9/32" orifice size, ASCO 8262152, for incoming effluent and dilution water pipes and mixing chamber pipe.
 - b. 1, normally closed, two-way, 15 psi, water, 3/8" pipe size, 3/8" orifice size, ASCO 8030865, for D-1 chamber evacuation pipe.
 - c. 12, normally closed, two-way, 36 psi, water, 1/4" pipe size, 9/32" orifice size. ASCO 8262C38, for remaining dilution chambers (D2-D6) and effluent chamber (E1-E6) evacuation pipes.
3. Stainless steel tubing, seamless, austenitic, 304 grade for freshwater and 316 grade for saline water.
 - a. 10 ft of 3/8" OD, 0.035" wall thickness, for dilution water and effluent pipes.
 - b. 60 ft of 1/4" OD, 0.035" wall thickness, for dilution water and effluent pipes.
 - c. 1 ft of 3/4" OD, 0.035" wall thickness, for standpipe in D1 chamber.
4. Swagelok tube connectors, stainless steel.
 - a. 4, male tube connectors, male pipe size 1/4", tube OD 3/8".
 - b. 2, male tube connectors, male pipe size 1/2", tube OD 3/8".
 - c. 26, male tube connectors, male pipe size 1/4", tube OD 1/4".
 - d. 2, male tube connectors, male pipe size 3/8", tube OD 3/8".
 - e. 2, male adaptor, tube to pipe, male size 1/2", tube OD 3/8".
5. 7, 1200 mL stainless steel beakers.
6. Several lbs each of Neoprene stoppers, sizes 00, 0, and 1; 1 lb of size 5.
7. 14 - aquarium (1-20 liters).
8. Magnetic stirrer.
9. 2 - PVC ball valves, 1/2" pipe size.
10. Dilutor control panel - see Fig. 9 and equipment list.
11. Plywood sheeting, exterior grade: one - 4' x 8' x 3/4", one - 4' x 8' x 1/2".
12. Pine or redwood board, 1" x 8", 20 ft.
13. Epoxy paint, 1 gal.
14. Assorted wood screws, nails, etc.
15. 25 ft - 14" ID, TEFLON[®] tubing, to connect the mixing chambers to the test chambers.




Figure 4. Vacuum siphon dilutor system, general diagram (not to scale).

VACUUM SIPHON SYSTEM EQUIPMENT LIST

1. Dilutor Glass.
2. Stainless steel solenoid valves.
 - a. 2, normally open, two-way, 55 psi, water, 1/4" pipe size, 9/32" orifice size, ASCO 8262152, for incoming effluent and dilution water pipes.
 - b. 2, normally closed, two-way, 15 psi, water, 3/8" pipe size, 3/8" orifice size, ASCO 8030865, for dilution water chamber D-6 and effluent chamber E-2.
3. Stainless steel tubing, seamless, austenitic, 304 grade for freshwater and 316 grade for saline water.
 - a. 60 ft of 3/8" OD, 0.035" wall thickness, for dilution water and effluent pipes.
 - b. 20 ft of 5/16" OD, 0.035" wall thickness, for standpipes in mixing chambers.
 - c. 1 ft of 3/4" OD, 0.035" wall thickness, for standpipe in D1 chamber.
4. Swagelok tube connectors, stainless steel.
 - a. 4, male tube connectors, male pipe size 1/4", tube OD 3/8".
 - b. 2, male tube connectors, male pipe size 3/8", tube OD 3/8".
 - c. 2, male adaptor, tube to pipe, male pipe size 1/2", tube OD 3/8".
 - d. 2, male tube connectors, male pipe size 1/2", tube OD 3/8".
5. 7, 1,200 mL stainless steel beakers.
6. Several lbs each of Neoprene stoppers, sizes 00, 0, and 1; 1 lb of size 5.
7. 14 - aquarium (1-20 liters).
8. Magnetic stirrer.
9. 2, PVC Ball valves, 1/2" pipe size.
10. Dilutor control panel equipment - see Fig. 9 and equipment list.
11. 7, 120 mL NALGENE® bottles.
12. 3 ft, 1-in-2 aluminum bar, for siphon support brackets.
13. Stainless steel set screws, box of 50, for securing SS tubing in siphon support brackets.
14. Stainless steel hose clamps, box of 10, size #4 or 5, (need 3 boxes).
15. 6, NALGENE® T's, 5/16" OD.
16. 12, TYGON® Y connectors, 3/8" I.D.
17. TYGON® tubing, 3/8" OD, 10 ft.
18. Plywood sheeting, exterior grade: one - 4' x 8' x 3/4", one - 4' x 8' x 1/2".
19. Pine or redwood board, 1" x 8", 20 ft.
20. Epoxy paint, 1 gal.
21. Assorted wood screws, nails, etc.
22. 25 ft of 5/16" ID, TEFLON® tubing, to connect the mixing chambers to the test chambers.





Figure 5. Vacuum siphon dilutor system, detailed diagram (not to scale).



DRAIN HOLES IN BOTTOM PLATE (C1 AND C2) SHOWN FOR SOLENOID VALVE DILUTOR SYSTEM. FOR VACUUM SIPHON DILUTOR SYSTEM, DRAIN HOLE IS REQUIRED ONLY FOR CHAMBER E1.

INDIVIDUAL PART SIZE AND NUMBER OF PIECES USING 6 mm (1/4 in.) PLATE GLASS. NOTE: 1.6 mm (1/16 in.) No. 304 GRADE (FOR FRESH WATER) OR No. 316 GRADE (FOR SALINE WATER) STAINLESS STEEL MAY BE SUBSTITUTED FOR GLASS.


LENGTH WIDTH NO. PIECES (16)

A:	251 mm x 95 mm	4 (END PLATES)
B:	208 mm x 95 mm	5 (PARTITIONS)
C1:	107 mm x 107 mm	1 (BOTTOM PLATE FOR E1)
C2:	447 mm x 107 mm	1 (BOTTOM PLATE FOR D1-D6)
F1:	107 mm x 251 mm	2 (FRONT AND BACK PANELS FOR E1)
F2:	447 mm x 251 mm	2 (FRONT AND BACK PANELS FOR D1-D6)

INSIDE CELL MEASUREMENTS AND APPROXIMATE VOLUMES

	WIDTH	LENGTH	HEIGHT	VOLUME
E1:	95 mm x 95 mm	x 251 mm	-	2365 mL
D1:	125 mm x 95 mm	x 200 mm	-	2375 mL
D2:	40 mm x 95 mm	x 200 mm	-	780 mL
D3:	50 mm x 95 mm	x 200 mm	-	980 mL
D4:	60 mm x 95 mm	x 200 mm	-	1140 mL
D5:	60 mm x 95 mm	x 200 mm	-	1140 mL
D6:	70 mm x 95 mm	x 200 mm	-	1330 mL

Figure 6. Effluent and dilution water chambers (not to scale).



DRAIN HOLES IN BOTTOM PLATE (C) SHOWN FOR SOLENOID VALVE DILUTOR SYSTEM ONLY. FOR VACUUM SIPHON DILUTOR SYSTEM, A DRAIN HOLE IS REQUIRED ONLY FOR CHAMBER E2.


INDIVIDUAL PART SIZE AND NUMBER OF PIECES USING A 6 mm (1/4 in.) PLATE GLASS ARE SHOWN BELOW. NOTE: 1/16 in. No. 304 (FOR FRESH WATER) OR No. 316 STAINLESS STEEL (FOR SALINE WATER) MAY BE SUBSTITUTED FOR GLASS.

	LENGTH	WIDTH	NO. PIECES (8)
A	180 mm x	80 mm	= 2 (END PLATES)
B	155 mm x	80 mm	= 4 (PARTITIONS)
C	296 mm x	92 mm	= 1 (BOTTOM PLATE)
D	296 mm x	180 mm	= 2 (FRONT AND BACK PLATES)

INSIDE CHAMBER MEASUREMENTS AND APPROXIMATE VOLUMES.


	WIDTH	LENGTH	HEIGHT	VOLUME
E2:	110 mm x	80 mm	x 155 mm	= 1364 mL
E3:	60 mm x	80 mm	x 155 mm	= 744 mL
E4:	30 mm x	80 mm	x 155 mm	= 372 mL
E5:	30 mm x	80 mm	x 155 mm	= 372 mL
E6:	30 mm x	80 mm	x 155 mm	= 372 mL

Figure 7. Effluent chambers (not to scale).



SIDE VIEW

END VIEW



INDIVIDUAL PART SIZE AND NUMBER OF PIECES USING
6 mm (1/4 in.) PLATE GLASS. APPROXIMATE CAPACITY
4360 mL

M-1	125	mm	x	153	mm	-	1	(END PLATE, WITH HOLE)
M-2	125	mm	x	153	mm	-	1	(END PLATE)
M-3	240	mm	x	165	mm	-	1	(BOTTOM PLATE)
M-4	240	mm	x	125	mm	-	2	(SIDE PLATES)

Figure 8. Pre-mixing chamber (not to scale).




Figure 9. Dilutor control panel wiring diagram.

DILUTOR CONTROL PANEL EQUIPMENT LIST*

Designation	CKT Description	Manufacturer
A ₁	Encapsulated amplifier	Cutler Hammer 13535H98C
CTR-1	Cycle counter	Redington #P2-1006
ET	Elapsed time indicator	Conrac #636W-AA H&T
F ₁	Input power fuse	Little fuse 342038
J ₁	Receptacle	Amphenol 91PC4F
J ₂	Aux A.C. output jack	Stand. 3-prong AC Rept.
J ₃	Main input power cord	Stand. 3-prong AC male plug
L ₁	Fill indicator light	Dialco 95-0408-09-141
L ₂	Emptying indicator light	Dialco 95-0408-09-141
L.S.	Liquid level sensor (Dual Sensing Probe)	Cutler Hammer 13653H2
P ₁	Plug	Amphenol 91MC4M
S ₁	On-off main power switch (spst)	Cutler Hammer 7580 K7
S ₂	On-off aux power switch (spst)	Cutler Hammer 7580 K7
SJ ₁	Solenoid	(See Solenoid and Vacuum System equipment lists)
SJ ₂	"	" " "
SJ ₃	"	" " "
SJ ₄ -SJ ₆	Additional Solenoids for Solenoid Valve System	" " "
TDR-1	Time delay relay	Dayton 5x829
TDR-2	Aux time delay relay	Dayton 5x829

*Consult local electric supply house.

APPENDIX D

PLANS FOR MOBILE TOXICITY TEST LABORATORY

D.1. TANDEM-AXLE TRAILER




Figure 1. Mobile bioassay laboratory, tandem axle trailer. Above - external side view; below - internal view from above.




Figure 2. Mobile bioassay laboratory, tandem-axle trailer, external and internal end views.




Figure 3. Mobile bioassay laboratory, tandem-axle trailer, internal views of side walls.

APPENDIX D

PLANS FOR MOBILE TOXICITY TEST LABORATORY

D.2. FIFTH WHEEL TRAILER




Figure 4. Mobile bioassay trailer, fifth-wheel trailer, internal view from above.

APPENDIX E

CHECK LISTS AND INFORMATION SHEETS

E.1. TOXICITY TEST FIELD EQUIPMENT LIST

Truck

☐ Boards
☐ Cinder blocks
☐ Drums: ☐ 500 gal nalgene
 ☐ 55 gal metal - diesel fuel
 ☐ 22 gal
☐ Gas can ☐ 5 gal
☐ Jacks
☐ Jumper cables
☐ Oil
☐ Pumps: ☐ (2) Homelite
 ☐ Hoses & couplings
☐ Shovels
☐ Spare tires (trailer, generator)

☐ Brine shrimp eggs
☐ Broom
☐ Brushes (wash)
☐ Buckets
☐ Camera
☐ Chlorine kit (w/chem)
☐ Cleanser
☐ Clip board (lg, sm)
☐ Cork borer set
☐ Culture dishes (200 mL, Daphnia)
☐ Daphnia food
☐ Data sheets: ☐ Bioassay (static)
 ☐ Bioassay (flow-thru)
 ☐ Dilutor volume delivery
 ☐ Calibrator delivery sheet
 ☐ Daily events log

Trailer

☐ Acetone
☐ Aerators (battery operated)
☐ Air line: ☐ Clamps
 ☐ Aerators (battery operated)
 ☐ Air line: ☐ Clamps
 ☐ Stones
 ☐ Tubing
 ☐ Valves

☐ Alcohol
☐ Aluminum foil
☐ Alkalinity analysis (0.02 N H₂SO₄)
☐ Boots: ☐ safety
 ☐ wading
☐ Batteries ☐ D cell
☐ Beakers: ☐ 150 mL nalgene
 ☐ 200 mL glass (3 boxes)
☐ Bottles: ☐ D.O.
 ☐ wash
 ☐ Sample
 ☐ VOA vials
 ☐ 500 mL plastic
 ☐ Glass organic
 ☐ Qt. w/teflon liner

☐ Dish pan
☐ Dish rack
☐ Dissolved oxygen:
 ☐ KCL membrane solution
 ☐ Membranes
 ☐ Meter (YSI)
 ☐ Probes
 ☐ Reagent: ☐ MnSO₄
 ☐ Alkaline azide
 ☐ H₂SO₄
 ☐ 0.0375 Na thiosulfate
 ☐ Starch
☐ Distilled H₂O
☐ Emergency road kit
☐ Enamel pans (lg, sm)
☐ Erlenmeyer flasks: ☐ 500 mL (2)
 ☐ 1000 mL
 ☐ 2000 mL
☐ Extension cords
☐ Fire extinguisher
☐ First aid kit
☐ Fish nets, (lg, sm)
☐ Flash light
☐ Generator: ☐ Oil
 ☐ Filter - fuel
 ☐ Funnel
 ☐ Grease gun (wheels)
 ☐ Credit card
 ☐ Lock/key
 ☐ Siphon hose

E.1. TOXICITY TEST FIELD EQUIPMENT LIST (CONT.)

- | | |
|--|---|
| <p> <input type="checkbox"/> Glass cutter
 <input type="checkbox"/> Gloves (plastic)
 <input type="checkbox"/> Graduated cylinders:
 25 mL, 50 mL, 100 mL
 250 mL, 500 mL, 1000 mL, 2000 mL
 <input type="checkbox"/> Ground wire & rod
 <input type="checkbox"/> Hand soap
 <input type="checkbox"/> Hard hats
 <input type="checkbox"/> Hardness analysis: <input type="checkbox"/> Buffer
 <input type="checkbox"/> EDTA
 <input type="checkbox"/> indicator

 <input type="checkbox"/> HCl (20%)
 <input type="checkbox"/> Heaters: <input type="checkbox"/> Aquarium
 <input type="checkbox"/> Space
 <input type="checkbox"/> Hose: <input type="checkbox"/> Clamps
 <input type="checkbox"/> Connectors
 <input type="checkbox"/> Ice chests

 <input type="checkbox"/> Jars: <input type="checkbox"/> 750 mL (4 boxes)
 <input type="checkbox"/> 3 gal (glass) (1)
 <input type="checkbox"/> 5 gal (glass) (1)
 <input type="checkbox"/> Sample jugs (2)
 <input type="checkbox"/> Kimwipes (lg, sm)
 <input type="checkbox"/> Lab coats (2)
 <input type="checkbox"/> Level
 <input type="checkbox"/> Light 110 V
 <input type="checkbox"/> Log book
 <input type="checkbox"/> Magnetic stirrers: <input type="checkbox"/> Lighted
 <input type="checkbox"/> Other

 <input type="checkbox"/> Mop

 <input type="checkbox"/> Paper towels
 <input type="checkbox"/> Parachute cord
 <input type="checkbox"/> Parafilm
 <input type="checkbox"/> Pencils, pens
 <input type="checkbox"/> pH: <input type="checkbox"/> Meters, Orion
 <input type="checkbox"/> Meters, corning
 <input type="checkbox"/> Buffers, 4,7,10
 <input type="checkbox"/> Probes (extras)
 <input type="checkbox"/> Pipets: <input type="checkbox"/> Bulbs
 <input type="checkbox"/> Eyedroppers
 <input type="checkbox"/> Volumetric (1 mL, 5 mL, 10 mL)
 <input type="checkbox"/> Plastic bags
 <input type="checkbox"/> Quality assurance - SPCP
 <input type="checkbox"/> Rain gear
 <input type="checkbox"/> Reconstituted hard water
 <input type="checkbox"/> Refractometer
 <input type="checkbox"/> Respirators (cartridges) </p> | <p> <input type="checkbox"/> Rubber bands
 <input type="checkbox"/> Ruler
 <input type="checkbox"/> Safety glasses
 <input type="checkbox"/> Safety manual
 <input type="checkbox"/> Sample labels
 <input type="checkbox"/> Scissors
 <input type="checkbox"/> Screen bioassay cups
 <input type="checkbox"/> Sea salts
 <input type="checkbox"/> Separatory funnels & racks
 <input type="checkbox"/> Silent giants
 <input type="checkbox"/> Silicon sealant
 <input type="checkbox"/> Solenoids (spare)
 <input type="checkbox"/> Stainless steel tubing pieces
 <input type="checkbox"/> Standard Methods Hand Book
 <input type="checkbox"/> Stirring bars
 <input type="checkbox"/> Stoppers (assorted)
 <input type="checkbox"/> Submersible pumps: <input type="checkbox"/> lg, sm.
 <input type="checkbox"/> screens

 <input type="checkbox"/> Super ice
 <input type="checkbox"/> Tablets (paper)
 <input type="checkbox"/> Tape: <input type="checkbox"/> Cellophane
 <input type="checkbox"/> Color coded
 <input type="checkbox"/> Electrician
 <input type="checkbox"/> Masking
 <input type="checkbox"/> Nylon
 <input type="checkbox"/> Thermometers: <input type="checkbox"/> Dial
 <input type="checkbox"/> Glass
 <input type="checkbox"/> Tools (lock/key)
 <input type="checkbox"/> Tygon tubing, 1/8", 1/4", 3/8"
 <input type="checkbox"/> Volumetric flasks (1000 mL, 2000 mL)
 <input type="checkbox"/> WD40
 <input type="checkbox"/> Weigh boats
 <input type="checkbox"/> Wire tags </p> |
|--|---|

APPENDIX E

CHECK LISTS AND INFORMATION SHEETS

E.2. INFORMATION CHECK LIST FOR ON-SITE INDUSTRIAL
OR MUNICIPAL WASTE TOXICITY TEST

1. PRE-TRIP INFORMATION

Facility Name: _____

Address: _____

Phone number: _____

Plant Representative(s): _____

Names, Titles, Addresses of Company Personnel:

A. To Receive Correspondence: _____

B. To Receive Carbons: _____

Date of Notification Letter: _____

State Making Notification and Arrangements: _____

Special Plant Safety/Security Requirements for EPA Personnel to Observe:

Local Accommodation Recommendations: _____

Directions to Plant: _____

Availability of Power Hookups (three 20-amp, 110-V Circuits): _____

Distance from Power Source to Trailer: _____ (Feet)

Trailer Location: _____

Possible Source of Dilution Water: _____

Major Products: _____

Raw Materials: _____

Name of Receiving Water: _____

Schedule of Plant Operation (continuous, weekdays only, etc.): _____

Treatment Steps: _____

Treatment Level (BPT, BAT, etc.): _____

Wastewater Retention Time by Lagoon or Treatment Step:

<u>Lagoon Designation</u>	<u>Retention Hours</u>	<u>Time Days</u>
_____	_____	_____
_____	_____	_____
_____	_____	_____
_____	_____	_____
_____	_____	_____

Total Wastewater Retention Time: _____ Hours; _____ Days

Retention Time Determination: _____ Calculated; _____ Actual

Calculation method: _____

Description of Wastewater Tap Point: _____

Description of Outfall (surface, submerged diffuser, etc.): _____

Description of Other Waste Disposal Alternatives in Use (spray irrigation, deepwell, municipal discharge, etc.): _____

2. ON-SITE INFORMATION

Wastewater General Characteristics:

Color: _____

Odor: _____

Solids: _____

Other: _____

Serial Number(s) of Discharge(s) to be Tested: _____

Description of Receiving Water: _____ Uniflow; _____ Tidal;
_____ Approximate amplitude, feet

Color: _____

Odor: _____

Solids: _____

Salinity: High tide _____; Low tide _____

Other: _____

7Q10: _____; Ave. flow _____

Description of Receiving Water Zone of Dilution: _____

Location and Description of Water Sampling Point(s): _____

Fresh: _____

Salt: _____

Dilution Waste General Characteristics:

Color: _____

Odor: _____

Solids: _____

Other: _____

Description of Toxicity Test Anomalies (plant production changes, power failure, rain events, etc.):

Duration		Anomaly
Time & Date	Time & Date	
_____	_____	_____
_____	_____	_____
_____	_____	_____
_____	_____	_____
_____	_____	_____
_____	_____	_____

Description of Plant maintenance: _____

Attach: DIAGRAM OF WASTEWATER TREATMENT FACILITIES.

3. FOLLOW-UP INFORMATION

Date of follow-up letter: _____

Wastewater Flow (data supplied by discharger):

<u>Week Prior to Testing</u>		<u>Week of Testing</u>	
<u>Date</u>	<u>Discharge (MGD)</u>	<u>Date</u>	<u>Discharge (MGD)</u>
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____

Average Discharge (MGD): _____

Organisms Tested On-site or In-Lab:

<u>Species</u>	<u>Flow-thru test duration (h)</u>	<u>Static test duration (h)</u>	<u>Test Location</u>	<u>Dates</u>	<u>Results</u>
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____

Possible Recommended Action as a Result of These Findings:

APPENDIX E

CHECK LISTS AND INFORMATION SHEETS

E.4. DILUTOR CALIBRATION FORM

Calibration Site: _____ Calibrator: _____
 Dilutor Number: _____ Date: _____

Effluent Concentration (%)	100.0	50.0	25.0	12.5	6.25	3.12	1.56
Dilution Water (mL)	0	500	750	876	938	969	984
Trial 1							
2							
3							
Average							
Effluent (mL)	1000	500	250	125	62	31	16
Trial 1							
2							
3							
Average							

Mixing Chamber (%): _____
 Wastewater (mL): _____
 Dilution Water (mL): _____

	Dilution Water	Effluent
Vol (mL)		
Trial 1		
2		
3		
Average		

Remarks:

